

ANALYSIS FOR ISOCYANATE CONTENT

1. SCOPE: This method is designed for the determination of the isocyanate content of the resin system found in polyurethane paints. The reported isocyanate content is based upon 100% resin solids.
2. APPARATUS AND MATERIALS:
 - 2.1. Analytical balance
 - 2.2. High speed centrifuge
 - 2.3. Centrifuge tubes, minimum 50 mL heavy wall glass with vinyl lined screw caps
 - 2.4. 10 ml Syringes with Luer slip tip and caps
 - 2.5. 60 mm Aluminum weighing dishes
 - 2.6. Drying oven $110 \pm 5^{\circ}\text{C}$
 - 2.7. Automatic titrator with stirrer (Mettler-Toledo DL 50 or equivalent)
 - 2.8. pH electrodes
 - 2.9. Fleaker or other vessel with stopper or cap, 200 mL minimum volume
 - 2.10. o-Xylene, reagent grade, dried over 3A Molecular Sieve (See 4.1)
 - 2.11. Buffer solutions, pH 4.00 and pH 7.00
 - 2.12. Dibutylamine, 99% assay
 - 2.13. Methanol, ACS reagent grade, dried over 3A Molecular Sieve (See 4.1)
 - 2.14. Hydrochloric acid solution, 1.0 N

3. PROCEDURE:

- 3.1. To analyze a one-component material, the resin system must first be separated from the whole paint.
- 3.2. Place a 40 ml sample of the whole paint into a high-speed centrifuge tube and cap immediately.
- 3.3. Centrifuge the sample to produce a minimum relative centrifugal force (RCF) of 4300 until the pigment and resin have separated. See 5.1 for RCF equation. See 4.2 for products requiring dilution to facilitate separation.
- 3.4. No separation is required for multi-component materials, when the resin system is packaged separately. The resin system is generally labeled “curing agent” or “part B”.
- 3.5. Determine the percent solids by weight of the resin solution (the resin solution is considered to be either the supernatant of a one-component material or the resin constituent of a multi-component material).
- 3.6. Condition the 60 mm aluminum weighing dishes in a drying oven at $110 \pm 5^{\circ}\text{C}$ for a minimum of 10 minutes and cool to ambient conditions in a desiccator.
- 3.7. Weigh 0.5 ± 0.1 g to 1.0 mg of the resin solution by difference from a syringe into a conditioned 60 mm aluminum weighing dish.
- 3.8. Immediately place the sample into a drying oven at $110 \pm 5^{\circ}\text{C}$ for 1.0 ± 0.25 hours.
- 3.9. Allow sample to cool to ambient conditions in a desiccator.
- 3.10. Weigh sample to 1.0 mg.
- 3.11. Calculate the percent solids by weight of the resin solution. See 5.2 for percent solids calculation.
- 3.12. Calibrate the titrimeter and pH electrodes with standard pH 7.00 and pH 4.00 solutions.
- 3.13. Place 25 ml of dry o-xylene into an appropriate size fleaker or vessel.
- 3.14. Weigh in the appropriate amount of resin solution. See 5.3 for sample weight equation.
- 3.15. Record the sample weight to the nearest 1.0 mg.
- 3.16. Add 20.00 ml of 1.7 N dibutylamine solution. Prepared by diluting 290 ml of dibutylamine to 1.0 L with dried o-xylene.

- 3.17. Cap the container immediately.
- 3.18. Allow the sample to stir for a minimum of 30 minutes but no longer than 60 minutes at room temperature.
- 3.19. Prepare a blank in the same manner, omitting the resin solution.
- 3.20. Wash down the sides of the containers with 100 mL of methanol.
- 3.21. Titrate the stirred blank with 1.0 N hydrochloric acid through the end point, which occurs at pH 5.0.
- 3.22. Record the volume of titrant used to the nearest 0.01 mL.
- 3.23. Titrate the stirred sample with 1.0 N hydrochloric acid through the end point, which occurs at pH 5.0.
- 3.24. Record the volume of titrant used to the nearest 0.01 mL.
- 3.25. Calculate the weight percent NCO. See 5.4 for weight percent NCO calculation.

4. NOTES:

- 4.1. Methanol and o-xylene used for this titration should be passed through separate columns containing 3A molecular sieve to remove any moisture that may be present. A recommended flow rate for passing methanol or o-xylene through a drying column is approximately 100 ml per minute. Dried stock may be collected into amber bottles and preserved by placing a small amount of 3A molecular sieve in the bottles.
- 4.2. Separation of pigment and resin of one-component materials may be aided by diluting the material up to the maximum volume of the centrifuge tube using dry o-xylene. Seal the centrifuge tube and shake the sample to obtain a homogenous sample prior to placing sample in the centrifuge.

5. CALCULATIONS:

- 5.1. Calculation for Relative Centrifugal Force (RCF):

$$RCF = 0.00001118 \times r \times N^2$$

Where, RCF = Relative Centrifugal Force

r = Rotating radius (cm)

N = Revolutions per minute

5.2. Calculation of weight percent resin solids of the resin solution:

$$S_R = (W_2 / W_1) \times 100$$

Where, S_R = Weight percent solids of the resin solution

W_2 = Weight of sample after drying

W_1 = Initial weight of sample

5.3. Calculation of the appropriate amount of resin sample to be used:

$$W = (42 / \text{NCO}) / R$$

Where, W = the weight in grams of the resin solution to be used

42 = the equivalent weight of an isocyanate group

NCO = the expected % NCO in the resin system

R = the percent weight solids of the resin solution

5.4. Calculation of the weight percent NCO of the sample:

$$\% \text{NCO} = (((B-S) \times N \times 42) / (1000 \times (W \times R))) \times 100$$

Where, B = the amount of titrant used to titrate the blank in milliliters

S = the amount of titrant used to titrate the sample in milliliters

N = the normality of the hydrochloric acid solution in equivalents per liter

42 = the equivalent weight of an NCO group

1000 = the factor converting liters to milliliters

W = the weight of the resin solution sample in grams

R = the percent weight solids of the resin solution

6. REPORT:

6.1. Report the percent weight solids of the resin solution to the nearest 0.01%.

6.2. Report the resin solution sample weight to the nearest milligram, 0.001 g.

6.3. Report the percent NCO of the resin system to the nearest 0.1%.

APPROVED _____
Director
DIVISION OF MATERIALS

DATE 1/7/03

Kentucky Method 64-250-03
Revised 1/7/03
Supersedes 64-250-2000
Dated 5/8/2000

Km25003.doc